



Identify Tooth and Shade in the Diagram to left

Ditch dies to identify margins. Do not over ditch.



Send pinned and sectioned modelwork with dies ditched similar to the illustration. Do not paint or treat with die hardener. Do not block-out with wax or any other glossy substance.

ITEMS ENCLOSED IN SHIPMENT

- Model Upper
Model Lower
Single Die
- Impression
Bite
Articulator

Please indicate Coping or Full Countour for each tooth

☐ Coping

☐ Full Contour

RESTORATION INSTRUCTIONS

Custom Abutment

Patient Name

Date Due

Additional Instructions

Lab Name

Contact

Address

Phone

Shipping Method:

2nd Day
Next Day Saver
Next Day Air